

Doppler-Free Spectroscopy

Rosalyn Chan

York University

PHYS4160

Cody Storry

16th November 2021

Introduction

The Doppler-free spectroscopy experiment demonstrates the general structures of the rubidium atom with reference to quantum numbers n , L , S , J and F . The origin of cross-over peaks and lamb dips are defined and described as well as determining the frequency difference of the rubidium hyperfine splitting of the $5P_{3/2}$ state.

Section 6.3.2: Pump beam alignment

Think about it:

Given the selection rules, the allowable transitions must satisfy $\Delta L = 1$ and $\Delta F = -1, 0, +1$. Therefore, the origin of Lamb dips we would expect to observe for the $^{85}\text{Rb } 5S_{1/2} F = 3$ to $5P_{3/2}$ states would be from $5S_{1/2} F = 3$ to $5P_{3/2} F = 2$, $5S_{1/2} F = 3$ to $5P_{3/2} F = 3$, $5S_{1/2} F = 3$ to $5P_{3/2} F = 4$ and also exactly halfway in between every pair of allowable transitions of $5S_{1/2} F = 3$ to $5P_{3/2} F = 2$, $5S_{1/2} F = 3$ to $5P_{3/2} F = 3$, $5S_{1/2} F = 3$ to $5P_{3/2} F = 4$.

Section 6.3.3: Background subtraction

Think about it:

The background subtracted beam needs to be pass through the Rubidium cell for this experiment as this would result in the background beam being doppler broadened as well. Hence, when the background subtraction is performed, this would very nicely produce a visual where we are able to clearly see the hyperfine splitting.

Data Analysis

Oscilloscope traces for four Doppler-broadened absorption features of both the laser passing through the cell and the background beam is recorded with range frequency at position 5 and different amplitude settings for each of the four Doppler-broadened absorption features.

(a) Part 1:

The rate at which the laser is scanning its frequency back and forth is of range frequency at position 5. The inverse of this rate is the amount of time the laser takes to undergo a full scan cycle. Half of this value will be the amount of time the laser takes to ramp up. In order to find the spacing of two peaks in units of laser frequency, the spacing between two peaks is determined in seconds and divide it by the number of seconds it takes for the full ramp up and then multiplying this ratio into the scan amplitude calibration.

Amount of time the laser takes to undergo a full scan cycle: $0.01\text{s} \pm 0.00005\text{s}$

Amount of time the laser takes to ramp up: $0.005\text{s} \pm 0.000025\text{s}$

The uncertainty of time is taken to be 0.00005s , that is by taking the lowest division of the oscilloscope scale.

^{87}Rb F=1 to F'=0,1,2:

Amplitude setting: 3.5

Amplitude calibration at setting 3.5: $6.514 \times 10^9 \text{Hz}$




Figure 1: Lamb dips of ^{87}Rb F=1 at amplitude setting 3.5.

^{85}Rb F=2 to F'=1,2,3:

Amplitude setting:1.5

Amplitude calibration at setting 1.5: $4.545 \times 10^9 \text{Hz}$




Figure 2: Lamb dips of ^{87}Rb F=2 at amplitude setting 1.5.

^{85}Rb F=3 to F'=2,3,4:

Amplitude setting:0.5

Amplitude calibration at setting 0.5: $1.970 \times 10^9 \text{Hz}$




Figure 2: Lamb dips of ^{85}Rb F=3 at amplitude setting 0.5.

^{87}Rb F=2 to F'=1,2,3:

Amplitude setting:1.0

Amplitude calibration at setting 1.0: $3.939 \times 10^9 \text{ Hz}$




Figure 4: Lamb dips of ^{87}Rb $F=2$ at amplitude setting 1.

Part 2:

The peaks are labelled by its corresponding excited state F number. The crossover peaks are labelled with a pair of numbers which denotes which two excited states contributes to their formation. For instance, a peak labelled as (2,3) is the crossover between $F'=2$ and the $F'=3$ excited states.

^{87}Rb F=1 to F'=0,1,2




Figure 5: Lamb dips of ^{87}Rb F=1 at amplitude setting 3.5.

No lamb dips can be seen in ^{87}Rb F=1 to F'=0,1,2

^{85}Rb F=2 to F'=1,2,3




Figure 6: Lamb dips of ^{85}Rb F=2 at amplitude setting 1.5.

The origin of the lamb dips would be from $5S_{1/2}$ F=2 to $5P_{3/2}$ F'=2.

^{85}Rb F=3 to F'=2,3,4




Figure 7: Lamb dips of ^{85}Rb F=3 at amplitude setting 0.5.

The origin of the lamb dips would be from $5S_{1/2}$ F=3 to $5P_{3/2}$ F'=3.

^{87}Rb F=2 to F'=1,2,3




Figure 8: Lamb dips of ^{87}Rb F=2 at amplitude setting 1.0

The origin of the lamb dips would be from $5S_{1/2} F=2$ to $5P_{3/2} F'=2$.

Part 3:

The hyperfine splitting of $5P_{3/2}$ F states is determined by computing the difference of two lamb dips in seconds, and dividing by the amount of time the laser took to ramp up and multiplying that ratio into the amplitude calibration.

The uncertainty of time is taken to be 0.00005s, that is by taking the lowest division of the oscilloscope scale.

^{87}Rb F=1 to F'=0,1,2

Amplitude setting: 3.5

Amplitude calibration: $6.514 \times 10^9 \text{Hz}$

No lamb dips can be seen in ^{87}Rb F=1 to F'=0,1,2, therefore the hyperfine splitting can not be carried out as one requires to identify the lamb dips correctly.

^{85}Rb F=2 to F'=1,2,3

Amplitude setting: 1.5

Amplitude calibration: $4.545 \times 10^9 \text{Hz}$

Time between two lamb dips (s), t_1 : $0.0138\text{s} - 0.00936\text{s} = 4.44 \times 10^{-3} \text{s} \pm 0.00005\text{s}$

Time of the laser going through a single ramp, t_2 : $0.005\text{s} \pm 0.000025\text{s}$

Ratio of time, t_r : $\frac{4.44 \times 10^{-3}}{0.005} = 0.888$

Uncertainty of the time ratio $\delta_t = \sqrt{\left(\frac{\delta_{t_1}}{t_1}\right)^2 + \left(\frac{\delta_{t_2}}{t_2}\right)^2} = \sqrt{\left(\frac{0.00005}{4.44 \times 10^{-3}}\right)^2 + \left(\frac{0.000025}{0.005}\right)^2} = 0.0123$

\Rightarrow time ratio, t_r : 0.888 ± 0.012

\therefore Hyperfine splitting, H: $\left(\frac{4.44 \times 10^{-3}}{0.005} \pm 0.012\right)(4.545 \times 10^9) \text{Hz} = (4.036 \pm 0.055) \times 10^9 \text{Hz}$

^{85}Rb F=3 to F'=2,3,4

Amplitude setting: 0.5

Amplitude calibration: $1.970 \times 10^9 \text{Hz}$

Time between two lamb dips, t_1 : $0.01672\text{s} - 0.01128\text{s} = 5.44 \times 10^{-3} \text{s} \pm 0.00005\text{s}$

Time of the laser going through a single ramp, t_2 : $0.005\text{s} \pm 0.000025\text{s}$

Ratio of time, t_r : $\frac{5.44 \times 10^{-3}}{0.005} = 1.088$

Uncertainty of the time ratio $\delta_t = \sqrt{\left(\frac{\delta_{t_1}}{t_1}\right)^2 + \left(\frac{\delta_{t_2}}{t_2}\right)^2} = \sqrt{\left(\frac{0.00005}{5.44 \times 10^{-3}}\right)^2 + \left(\frac{0.000025}{0.005}\right)^2} = 0.011$

\Rightarrow time ratio, t_r : 1.088 ± 0.011

\therefore Hyperfine splitting: $\left(\frac{5.44 \times 10^{-3}}{0.005} \pm 0.011\right)(1.970 \times 10^9) \text{Hz} = (2.143 \pm 0.022) \times 10^9 \text{Hz}$

^{87}Rb F=2 to F'=1,2,3

Amplitude setting: 1.0

Amplitude calibration: $3.939 \times 10^9 \text{Hz}$

Time between two lamb dips, t_1 : $0.0156\text{s} - 0.01244\text{s} = 3.16 \times 10^{-3} \text{s}$

Time of a laser going through a single ramp, t_2 : 0.005s

Ratio of time, t_r : $\frac{3.16 \times 10^{-3}}{0.005} = 0.632$

Uncertainty of the time ratio $\delta_t = \sqrt{\left(\frac{\delta_{t_1}}{t_1}\right)^2 + \left(\frac{\delta_{t_2}}{t_2}\right)^2} = \sqrt{\left(\frac{0.00005}{3.16 \times 10^{-3}}\right)^2 + \left(\frac{0.000025}{0.005}\right)^2} = 0.016$




\Rightarrow time ratio, t_r : 0.632 ± 0.012

$$\therefore \text{Hyperfine splitting: } \left(\frac{3.16 \times 10^{-3}}{0.005} \pm 0.012 \right) (3.939 \times 10^9) \text{Hz} = (2.489 \pm 0.047) \times 10^9 \text{Hz}$$


Discussion and Conclusion:

The hyperfine splitting of the $5P_{3/2}$ F states are found to be $(4.036 \pm 0.055) \times 10^9 \text{Hz}$ for ^{85}Rb F=2, $(2.143 \pm 0.022) \times 10^9 \text{Hz}$ for ^{85}Rb F=3 and $(2.489 \pm 0.047) \times 10^9 \text{Hz}$ for ^{87}Rb F=2. The hyperfine splitting of ^{87}Rb F=1 could not be computed as no lamb dips were able to be identified. It is also important to note here that the hyperfine splitting found here are at an order of GHz which is much higher than the expected value at an order of Mhz. One of the many possibilities is due to the fact that not all six peaks were able to be produced during the whole process of the experiment. As a matter of fact, technical issues and complications were also ran upon during the duration of the experiment which might have contributed to the data being somewhat unreliable.

Appendix A:

Rubidium states	Oscilloscope traces	Amplitude setting
^{87}Rb F=1 to F'=0,1,2		3.5
^{85}Rb F=2 to F'=1,2,3		1.5
^{85}Rb F=3 to F'=2,3,4		0.5

^{87}Rb F=2 to F'=1,2,3



1.0